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

TexiShield® – Item reference: F2-B IIR

Instructions for use of respiratory masks FFP2 NR D / Type IIR for single use

1. Description and indication

Respiratory masks FFP2 NR D / Type IIR for single use from MDB Texinov are filtering half-mask to protect against particles, they meet the requirements of the standards:

- EN 149:2001+A1:2009: the masks are personal protective equipments (PPE) that protect against liquid or solid sprays (dust, smokes, aqueous aerosols...), usable by Industry professionals. They do not protect against gas and vapour.
- EN 14683+AC:2019 / ISO 22609:2004: the masks are Class I medical device (MD), usable by health professionals and patients, that protect against splashes of infectious agents.

Composition: polypropylene, spunbond and meltblown, nose clip with metal wires, elastomer fasteners.

2. Control before use

This mask is suitable for people trained for its use. Check the storage expiration date on the packaging. Don't use any mask with an expired storage date. Take the mask out of its pack at the time of using it only.

Before using it, wash your hands and ensure that it is intact and that the elasticised strings, that split in two parts, are intact too. If the mask is damaged, put it throw.

To ensure that the respiratory mask tightly adheres to your face, cup your hands over the mask to block up its surface of filtration as much as possible without moving it and breathe out strongly. If any air flows around your nose, readjust the nose bar. If any air flows around the mask, carefully adjust the strings to the back of the head.

When you breathe in, the mask should deform.

In case tightness obtained is not satisfactory, do not enter the area where it is required to wear the mask. Choose a different model of mask to obtain the required tightness.

3. Use

Follow instructions 1→2→3→4→5→6→7→8

This mask does not require any maintenance. It has been designed for one single continuous work session i.e. 8 hours. Any failure to comply with use and adjustment controls and warnings concerning the use of the respiratory mask for the whole period of exposure, may reduce its efficiency.

This mask is for single use.

The reuse of a mask already used may cause risks of contamination of the vector and people nearby because of the deterioration of the performances of the mask.

It's unlikely that the requirements for tightness are met if beard's hair go through the face seal.

Leave the premises and dispose the mask if:

- You feel dizzy or feel any pain
- It becomes hard to breathe: oxygene is missing
- Air quality is compromised (contaminants)
- It is damaged
- It seems to hamper the breathe
- It is tainted by blood or any other infectious substance.

Dispose of the mask according to the current regulation.

Do not use in explosive atmosphere.

Any serious incident in connection with the device should be notified to the manufacturer and to the competent authority of the Member State in which the user and / or patient is established.

4. Storage and Expiry date

Storage according to the recommendations on the packaging (between 5 to 40°C and below 85% of humidity).

The product expires 5 years after its date of manufacture. Expiry date is indicated on the mask and on the packaging.

5. Performance justifying the protection class

This filtrating half-mask meets the protection requirements defined by the regulation PPE 2016/425, EN 149:2001+A1:2009 and the Regulation (EU) 2017/745.

FFP2 NR D / Type IIR Regulation (EU) PPE 2016/425 Regulation (EU) 2017/745		
EN 149 : 2001 + A1 : 2009	Paraffine oil initial penetration	≤ 6 %
	NaCl initial penetration	≤ 6 %
	Breathe resistance at 30L/min	≤ 0,7 mbar
	Breathe resistance at 95 L/min	≤ 2,4 mbar
	Breathe resistance at 160 L/min	≤ 3 mbar
EN 14683 + AC : 2019 ISO 22609 : 2004	Bacterial filtration efficiency	≥ 98%
	Differential pressure	< 60 Pa/cm ²
	Splash resistance pressure	≥16 kPa
	Microbiological cleanliness	≤ 30 CFU/g

6. Details of the certification laboratory

Notified body implicated into the EU Type Examination (Module B and D) according to regulation PPE 2016/425/EU: Apave Exploitation France SAS (n°0082) - 6 Rue du Général Audran - 92412 COURBEVOIE cedex - FRANCE

The internet address where the EU declaration of conformity can be assessed is: www.texinov.com. **CE Mark 0082**

7. Meaning of symbols

	Manufacturer		Medical Device
	Protect from light and sunlight		Maximum relative moisture of storage conditions <85%
	End of storage date		Temperature range of storage conditions
	Batch number		Product reference
	For single use		Read instructions for use
	Resistant to Dolomite clogging		Non reusable
	Recyclable		Latex free Graphene free Phthalate free Bisphenol free

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